FORM FIT®

HYDROGEL CANALICULAR PLUG

DESCRIPTION

The OASIS Medical Form Fit® canalicular plug is designed to be inserted through the punctal opening in order to block tear drainage through the canaliculus. The Form Fit® canalicular plug is provided as a semi-rigid rod of polyvinyl pyrrolidone (PVP) based hydrogel measuring approximately 0.3mm diameter by 2.5mm long. Upon contact with fluids, the plug will slowly swell to approximately 3 times its initial size to form a soft, mucus-like gel.

Each Form Fit® canalicular plug is supplied sterile, in dry form, in a disposable inserter. The product is gamma sterilized and sterility is guaranteed unless the packaging is opened or otherwise compromised. The Expiration Date is provided on the label and the Form Fit® canalicular plug should not be used after this date.

INDICATIONS FOR USE

The OASIS Medical Form Fit® canalicular plug is intended for use in patients experiencing

dry eye symptoms such as redness, burning, reflex tearing, itching or foreign body sensations which can be relieved by occluding the punctum. The Form Fit® canalicular plug may be used in the treatment of dry eye syndrome and the dry eye components of various ocular surface diseases. When indicated, the Form Fit® canalicular plug may be used after ocular surgery to prevent complications due to dry eye and to enhance the retention of topical ocular medications on the eye. Patients experiencing dry eye related contact lens problems may also be aided by the Form Fit® canalicular plug.

CONTRAINDICATIONS

The OASIS Medical Form Fit® canalicular plug is contraindicated for patients experiencing epiphora, inflammation of the eye lid, and tearing secondary to dacryocystitis with mucopurulent discharge. If the patient experiences irritation, infection, or epiphora after insertion of the Form Fit® canalicular plug, the plug should be removed.

INSERTION

Anesthetize the area of the punctum with either a topical anesthetic instilled into the conjunctival sac or an anesthetic soaked, cotton tipped applicator held against the punctum for 30 seconds. A small amount of lubricant may be used on the punctum in order to facilitate insertion across the punctum.

Caution: Patency of the canaliculus should be determined prior to insertion of the Form Fit® plug.

The OASIS Medical Form Fit® canalicular plug comes preloaded in a disposable inserter. Verify that the white Form Fit® plug has not fallen out of the inserter during shipping. Place the brown tip portion of the inserter through the puncta and squeeze the inserter to push the plug from the inserter tip into the canaliculus. The top edge of the plug should lie below the punctum opening. Forceps may be used to further push the Form Fit® canalicular plug to obtain the desired placement.

REMOVAL

Should it be necessary to remove the Form Fit® canalicular plug, flush the mucus-like, gel plug from the canaliculus using a syringe filled with saline. Place an injection cannula on the tip of the syringe, position the tip of the injection cannula through the puncta, and gently flush the canaliculus with saline.

PRECAUTIONS

Do not resterilize the Form Fit® canalicular plug by any method. Use only sterile intraocular irrigation solutions when required. Do not store Form Fit® canalicular plugs in excessively humid environments.

During implantation avoid getting the plug wet as this will cause it to expand. Avoid touching the plug or inserter tip to the patient's eyelid prematurely.

Should an infection occur in the canaliculus, flush the mucus-like plug from the canaliculus and allow the infection to clear. The new Form Fit® plug may be inserted into the canaliculus after the infection clears.

Care should be used not to perforate the punctum or canaliculus during the insertion of the Form Fit® canalicular plug. Perforation may cause pain, increase the risk of infection, and cause the plug to fall out more easily due to an enlarged punctal opening or infection. If perforation occurs, delay the placement of the plug until the wound heals.

The Form Fit® canalicular plug may enhance the effects of ocular medications used on the eye. Depending on the type of medications being used, dosage may need to be reduced accordingly.

DO NOT REUSE. This is a single use, disposable, device that cannot be properly cleaned and re-sterilized.

Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

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