

OASIS®





***CE Marked**



CORNEAL COMFORT & PROTECTION

Provides ocular surface protection following surgery, injection and traumatic or non-traumatic corneal conditions.

Available in Quick Shield, 12 hour, 24 hour, and 72 hour dissolution times.

Diameter: 14.05mm

Compound base curve: 9mm

PRODUCT FEATURES

- Comfort of a thin and pliable collagen barrier
- Collagen certified suitable for use by EDQM
- Hydrates to 70% water content
- Protects the surface of the eye
- Lubricates with collagen while it gradually dissolves
- CE approved

Product Name	Duration	Ref #
SOFT SHIELD Quick Shield Collagen Shield 10pk	Less than 12 hours	7002
SOFT SHIELD Collagen Shield 10pk	12 hours non-patched, 24 hours patched	7012
SOFT SHIELD Collagen Shield 10pk	24 hours non-patched, 72 hours patched	7024
SOFT SHIELD Collagen Shield 10pk	72 hours non-patched	7072



***CE Marked**

Product Name	Description	Qty	Ref #
Surgical Marking Pen	Gentian Violet pre-inked pen	10 per box	0442





*CE Marked

PRODUCT FEATURES

- 2ml filled syringe with cannula
- HPMC (hydroxypropyl methylcellulose)
- Stores at room temperature



*CE Marked

PRODUCT FEATURES

- 1.5ml filled syringe with cannula
- HPMC (hydroxypropyl methylcellulose)
- Stores at room temperature



PRODUCT FEATURES

- 2ml filled syringe with cannula
- Hypromellose
- Stores at room temperature

To Order Call (844) 820 – 8940

VISCO SHIELD[®] HPMC VISCOELASTIC

VISCO SHIELD[®] HPMC Viscoelastic is a dispersive HPMC viscoelastic composed of a sterile, high molecular weight, hydroxypropyl methylcellulose solution for use in anterior segment surgical procedures including cataract extraction and intraocular lens implantation.

FOR EXPORT ONLY

QTY	Description	Ref #
6 per box	20,000 cps	5122
6 per box	8,000 cps	51082

VISCO SHIELD[®] *Topical*

VISCO SHIELD[®] Topical is preservative-free and designed to keep the cornea hydrated. It adheres to the cornea and forms a thin, optically neutral clear film.

FOR EXPORT ONLY

QTY	Description	Ref #
6 per box	8,000 cps	5200

VISCO SHIELD[®] *Topical Drops*

OPHTHALMIC LUBRICANT DROPS

VISCO SHIELD[®] Topical Drops is preservative-free and designed to keep the ocular surface lubricated and hydrated to prevent drying of the eye.

US DOMESTIC ONLY

QTY	Description	Ref #
6 per box	5,000 cps	TD5200-0



SURGICAL CUSTOMIZED SOLUTIONS

SURGEON'S CHOICE CARD

Create a custom surgical kit to efficiently meet your supply needs. Select from a variety of safety scalpels, cystotomes, cannulas, and more from the OASIS® product catalog.



- Customize by surgeon preference
- 2 options available
- Improve operating room efficiency
- Consistent procedural care
- Ease of purchasing from a single supply source

Contact your local OASIS Sales Consultant to request a quote or email:

customsurgicalkit@oasismedical.com

For more information, contact Customer Service at (844) 820-8940.


OPTION A



SURGICAL CUSTOMIZED SOLUTIONS

SURGEON'S CHOICE CANNULA CARD

OPTION B



BULK PACKAGING

- Economical – reduced price per unit
- Saves time – no need to open boxes
- Less storage space needed
- Less waste disposal
- Provided sterile
- Automatic monthly shipments available

**Available for cannulas, needles, cystotomes
& Premier Edge scalpels.**



Surgical disposable products are now offered in this cost-effective and convenient Customized Inventory Management Program. The product is provided bulk and sterile in plastic bags, to eliminate waste and save time. Automatic monthly shipments are available with your choice of delivery on the 1st or 15th of each month.

Contact your local OASIS® Sales Consultant to request a quote or email:

customsurgicalkit@oasismedical.com



SCALPELS & KNIVES

OASIS® carries a wide variety of high quality ophthalmic microsurgical knives from Feather and Premier Edge for use in cataract, refractive, retinal and glaucoma surgery.



Feather redefines sharp and combines it with elegantly designed precision in their broad line of scalpels. Feather reusable scalpels are designed for repeated use, can be resterilized at your surgical facility in compliance with CMS regulations, and are easily disposed of after multiple uses. Feather also offers both single use disposable and safety scalpels.

Reusable	R	5 per box	Sure-grip aluminum handle	Stainless steel micro blades	Sterile individual autoclavable cases	
Disposable	D	5 per box	Disposable plastic handle	Stainless steel micro blades	Sterile blister pack	
Safeshield	S	5 per box	Disposable plastic handle	Stainless steel micro blades	Sterile blister pack	Retractable safety guard



***CE Marked**

Premier Edge Microsurgical Knives provide trustworthy performance during surgical procedures with enhanced precision on the cutting tip and ergonomic handle design. Safety scalpels are designed for one-handed use.

Disposable	D	6 per box	Disposable plastic handle	Stainless steel micro blades	Sterile blister pack	
Safety Line	S	6 per box	Disposable plastic handle	Stainless steel micro blades	Sterile blister pack	Retractable safety guard



GROOVE



R

FEATHER CONTROLLED DEPTH GROOVE

- Designed for controlled depth incisions
- Sterile, 5 per box

Product Name: FEATHER® MICRO FEATHER Ophthalmic Scalpel with Aluminum Handle

Depth	Reusable
300 microns (0.3mm)	07330G ●
350 microns (0.35mm)	07335G ●
400 microns (0.4mm)	07340G ●
550 microns (0.55mm)	07355G ●
600 microns (0.6mm)	07360G ●



R

FEATHER GLAUCOMA GROOVE

- Designed for fixed cutting depth
- Sterile, 5 per box

Product Name: FEATHER® MICRO FEATHER Ophthalmic Scalpel with Aluminum Handle

Depth	Thickness	Reusable
400 microns (0.4mm)	0.1mm	06340G ●

IMPLANT



R

FEATHER WOUND ENLARGEMENT IMPLANT

- Sterile, 5 per box

Product Name: FEATHER® MICRO FEATHER Ophthalmic Scalpel with Aluminum Handle

Width	Bevel	Blade	Reusable
3.5mm	Bevel Up	Angled	07635BR ●
5.1mm	Bevel Up	Angled	07651BR ●



D

PREMIER EDGE IMPLANT

- Sterile, 6 per box

Width	Bevel	Blade	Disposable
4.0mm	Bevel Up	Angled	PE3440 ●



SLIT KERATOME



R



FEATHER SLIT KNIFE RIB & R

- Sterile, 5 per box

Product Name: FEATHER® MICRO FEATHER Ophthalmic Scalpel with Aluminum Handle

Width	Bevel	Reusable	
2.2mm	Bevel Up	0622RC	●
2.4mm	Bevel Up	0624RC	●
2.75mm	Bevel Up	06275RC	●
2.8mm	Bevel Up	0628RC	●
3.0mm	Bevel Up	0630RC	●

- Improved Operability: Length from blade tip to bend is shortened from 8mm (conventional model) to 6mm.
- Arrowhead-shaped curved edge reduces incision resistance.
- Both sides of knife are thinly grinded (RIB shape).



R

D

FEATHER ULTRA-SHARP SLIT/KERATOME

- Sterile, 5 per box

Product Name: FEATHER® MICRO FEATHER Ophthalmic Scalpel with Aluminum Handle


Product Name: FEATHER® MICRO FEATHER Disposable Ophthalmic Scalpel

Width	Bevel	Blade	Reusable	Disposable
1.4mm	Bevel Up	Angled	0614BUC ●	P-0614BUC ●
1.6mm	Bevel Up	Angled	0616BUC ●	P-0616BUC ●
1.8mm	Bevel Up	Angled	0618BUC ●	P-0618BUC ●
2.2mm	Bevel Up	Angled	0622BUC ●	P-0622BUC ●
2.4mm	Bevel Up	Angled	0624BUC ●	P-0624BUC ●
2.5mm	Bevel Up	Angled	0625BUC ●	
2.6mm	Bevel Up	Angled	06265BUC ●	
2.75mm	Bevel Up	Angled	06275BUC ●	P-06275BUC ●
2.8mm	Bevel Up	Angled	0628BUC ●	P-0628BUC ●
3.0mm	Bevel Up	Angled	0630BUC ●	P-0630BUC ●
3.2mm	Bevel Up	Angled	0632BUC ●	P-0632BUC ●





R **D** **S**



FEATHER SLIT

- 60° keratome angle. 45° angle. 0.25mm blade thickness. • Sterile, 5 per box

Product Name: FEATHER® MICRO FEATHER Ophthalmic Scalpel with Aluminum Handle

Product Name: FEATHER® MICRO FEATHER Disposable Ophthalmic Scalpel

Product Name: FEATHER® MICRO FEATHER SAFESHIELD Disposable Ophthalmic Scalpel

Width	Bevel	Blade	Reusable	Disposable	Safety
1.2mm	Bevel Up	Angled	7612BU ●		
2.5mm	Double	Angled	07625B ●		
2.65mm	Double	Angled	07627B ●		
2.8mm	Double	Angled	07628B ●	P-7628B ●	PSS-7628B ●
3.0mm	Double	Angled	07630B ●	P-7630B ●	
3.0mm	Bevel Up	Angled	07630BU ●		
3.2mm	Double	Angled	07632B ●	P-7632B ●	
3.2mm	Bevel Up	Angled	07632BU ●		



D

PREMIER EDGE SLIT

- Sterile, 6 per box

Width	Bevel	Blade	Disposable
1.8mm	Bevel Up	Angled	PE3818 ●
2.2mm	Bevel Up	Angled	PE3822 ●
2.4mm	Bevel Up	Angled	PE3824 ●
2.5mm	Bevel Up	Angled	PE3825 ●
2.65mm	Bevel Up	Angled	PE3826 ●
2.75mm	Bevel Up	Angled	PE3827 ●
2.85mm	Bevel Up	Angled	PE3828 ●
3.0mm	Bevel Up	Angled	PE3830 ●
3.2mm	Bevel Up	Angled	PE3832 ●
3.2mm	Bevel Up	Straight	PE3832S ●



SLIT KERATOME



D **S**


PREMIER EDGE THIN SLIT

- Sterile, 6 per box

Width	Bevel	Blade	Disposable	Safety
2.0mm	Bevel Up	Angled	PE4820-TL ●	PE4820-SL ●
2.2mm	Bevel Up	Angled	PE4822-TL ●	PE4822-SL ●
2.4mm	Bevel Up	Angled	PE4824-TL ●	PE4824-SL ●
2.65mm	Bevel Up	Angled	PE4826-TL ●	PE4826-SL ●
2.75mm	Bevel Up	Angled	PE4827-TL ●	PE4827-SL ●
2.85mm	Bevel Up	Angled	PE4828-TL ●	PE4828-SL ●
3.0mm	Bevel Up	Angled	PE4830-TL ●	PE4830-SL ●
3.2mm	Bevel Up	Angled	PE4832-TL ●	PE4832-SL ●



R **D**



FEATHER CLEAR CORNEA SLIT

- 60° keratome angle. 45° angle. Depth markings at 1.5mm & 2.0mm
- Sterile, 5 per box

Product Name: FEATHER® MICRO FEATHER Ophthalmic Scalpel with Aluminum Handle

Product Name: FEATHER® MICRO FEATHER Disposable Ophthalmic Scalpel

Width	Bevel	Thickness	Reusable	Disposable
2.8mm	Double	0.15mm	09628B ●	P-9628B ●
3.0mm	Double	0.15mm	09630B ●	P-9630B ●
3.2mm	Double	0.15mm	09632B ●	



D **S**

PREMIER EDGE CLEAR CORNEA SLIT

- Sterile, 6 per box

Width	Bevel	Blade	Disposable	Safety
2.2mm	Double	Angled	PE3122 ●	
2.5mm	Double	Angled	PE3125 ●	
2.65mm	Double	Angled	PE3126 ●	PE3126-SL ●
2.75mm	Double	Angled	PE3127 ●	PE3127-SL ●
2.85mm	Double	Angled	PE3128 ●	PE3128-SL ●
3.2mm	Double	Angled	PE3132 ●	



R

FEATHER PARALLEL SLIT

- 50° keratome angle. 45° blade angle. Chamfered side edges
- Sterile, 5 per box

Product Name: FEATHER® MICRO FEATHER Ophthalmic Scalpel with Aluminum Handle

Width	Bevel	Thickness	Reusable
2.2mm	Bevel Up	0.15mm	07622BUPS ●
2.4mm	Bevel Up	0.15mm	07624BUPS ●
2.5mm	Bevel Up	0.15mm	07625BUPS ●



D

PREMIER EDGE CLEAR CORNEA PARALLEL SLIT

- Sterile, 6 per box

Width	Bevel	Blade	Disposable
2.65mm	Double	Angled	PE3226 ●
2.75mm	Double	Angled	PE3227 ●
2.85mm	Double	Angled	PE3228 ●
3.0mm	Double	Angled	PE3230 ●




D

PREMIER EDGE CLEAR CORNEA TRAPEZOID SLIT

- Sterile, 6 per box

Width Mark	Bevel	Blade	Disposable
1.5mm - 1.7mm	Bevel Up	Angled	PE3303 ●
1.8mm - 2.2mm	Bevel Up	Angled	PE3301 ●
2.5mm - 3.0mm	Bevel Up	Angled	PE3302 ●

MINI BLADES



D

PREMIER EDGE MINI BLADES

- Sterile

Description	Disposable
Multi-sided, sharp all around, 6 per box	PE5700 ●
Rounded Tip, sharp on one side, 20 per box	PE6400 ●
Curved Tip, sharp on one side, 20 per box	PE6700 ●
Rounded Tip, sharp all around, 6 per box	PE6900 ●



CRESCENT TUNNEL



R **D**

FEATHER CRESCENT TUNNEL

- 45° angle
- Sterile, 5 per box

Product Name: FEATHER® MICRO FEATHER Ophthalmic Scalpel with Aluminum Handle

Product Name: FEATHER® MICRO FEATHER Disposable Ophthalmic Scalpel

Width	Bevel	Blade	Reusable	Disposable
2.2mm	Double	Angled	07922B ●	
2.6mm	Bevel Up	Angled	07926BU ●	P-7926BU ●
2.6mm	Bevel Down	Angled	07926BD ●	
2.6mm	Double	Angled	07926B ●	P-7926B ●



D **S**

PREMIER EDGE CRESCENT TUNNEL

- Sterile, 6 per box

Width	Bevel	Blade	Disposable	Safety
2.0mm	Bevel Up	Angled	PE3720A ●	PE3720A-SL ●
2.0mm	Bevel Up	Straight	PE3720 ●	

SCLERAL



R

FEATHER MULTISIDED SCLERAL

- Sterile, 5 per box

Product Name: FEATHER® MICRO FEATHER Ophthalmic Scalpel with Aluminum Handle

Bevel	Thickness	Reusable
Double	0.38mm	0780 ●

ROUND TUNNEL



R **D** **S**

FEATHER ROUND TUNNEL

- Sterile, 5 per box
- 0.25mm blade thickness

Product Name: FEATHER® MICRO FEATHER Ophthalmic Scalpel with Aluminum Handle

Product Name: FEATHER® MICRO FEATHER Disposable Ophthalmic Scalpel

Product Name: FEATHER® MICRO FEATHER SAFESHIELD Disposable Ophthalmic Scalpel

Width	Bevel	Blade	Thickness	Reusable	Disposable	Safety
2.2mm	Double	Angled	0.25mm	07722B ●		
2.6mm	Double	Angled	0.25mm	07726B ●	P-7726B ●	
3.0mm	Double	Angled	0.25mm	07730B ●		PSS-7730B ●



D

PREMIER EDGE ROUND TUNNEL

- Sterile, 6 per box

Width	Bevel	Blade	Disposable
2.2mm	Double	Angled	PE3622 ●
3.0mm	Double	Angled	PE3630 ●

"The delivery of first-rate surgical outcomes is of paramount importance to our team. We constantly strive to produce those outcomes in the most efficient manner possible. Cost containment in our supply chain allows us to conserve our resources and to redirect them where they matter most. This is critical for any ambulatory surgery center to thrive. We use FEATHER Reusable Scalpels as one means to achieve these goals. The product is superiorly made and this allows us, along with the proper care, to get the most use from each blade."

Anonymous Independent ASC

"The Iris Expander's stable, atraumatic quality is preferable to rigid alternatives because it's easier to use and doesn't require anterior chamber injection of triamcinolone acetonide (Triesence, Alcon). I can insert it in almost any eye, even an eye with a ruptured anterior capsule. The device is also cost effective. I never hesitate to use it to reduce the potential for surgical complications that could be vision threatening as well as costly."

George R. Reiss, MD, Eye Physicians & Surgeons of Arizona

*See the OASIS® Iris Expander on page 28.



STAB INCISION



R **D** **S**

FEATHER STAB INCISION

- Sterile, 5 per box

Product Name: FEATHER® MICRO FEATHER Ophthalmic Scalpel with Aluminum Handle

Product Name: FEATHER® MICRO FEATHER Disposable Ophthalmic Scalpel

Product Name: FEATHER® MICRO FEATHER SAFESHIELD Disposable Ophthalmic Scalpel

Angle	Blade	Reusable	Disposable	Safety
15°	Straight	0715 ●	P-715 ●	PSS-715 ●
22.5°	Straight	0722 ●	P-722 ●	
30°	Straight	0730 ●	P-730 ●	PSS-730 ●
45°	Straight	0745 ●	P-745 ●	



D **S**

PREMIER EDGE STAB INCISION

- Sterile, 6 per box

Angle	Bevel	Blade	Disposable	Safety
15°	Double	Straight	PE3015 ●	PE3015-SL ●
22.5°	Double	Straight	PE3022 ●	
30°	Double	Straight	PE3030 ●	
45°	Double	Straight	PE3045 ●	



D

PREMIER EDGE RESTRICTED DEPTH STAB INCISION

- Sterile, 6 per box

Angle	Depth	Disposable
15°	3.0mm	PE3015-3 ●
15°	5.0mm	PE3015-5 ●

MVR




R

FEATHER MVR

- Sterile, 5 per box.

Product Name: FEATHER® MICRO FEATHER Ophthalmic Scalpel with Aluminum Handle

Gauge	Bevel	Blade	Reusable
20g	Double	Straight	07520 ●
23g	Double	Straight	07523 ●



D S

PREMIER EDGE MVR

- Sterile, 6 per box.

Gauge	Bevel	Blade	Disposable	Safety
19g	Double	Straight	PE3519 ●	
20g	Double	Straight	PE3520 ●	PE3520-SL ●
23g	Double	Straight	PE3523 ●	
25g	Double	Straight	PE3525 ●	
19g	Double	Angled	PE3519A ●	
20g	Double	Angled	PE3520A ●	
23g	Double	Angled	PE3523A ●	

SIDE PORT



R

FEATHER PARACENTESIS

- Sterile, 5 per box.

Product Name: FEATHER® MICRO FEATHER SAFESHIELD Disposable Ophthalmic Scalpel

Width	Bevel	Blade	Reusable
1.0mm	Bevel Up	Straight	07610 ●
1.0mm	Bevel Up	Angled	07610BU ●



D

PREMIER EDGE SIDE PORT

- Sterile, 6 per box.

Width	Bevel	Blade	Disposable
1.0mm	Double	Angled	PE3010 ●
1.2mm	Double	Angled	PE3012 ●
1.4mm	Double	Angled	PE3014 ●




DISPOSABLE MICROSURGICAL INSTRUMENTS

OASIS® supplies an extensive line of cannulas, needles, and cystotomes.

*CE Marked


ANESTHESIA



RETROBULBAR ANESTHESIA NEEDLE

- Sterile, 10 per box


Description	Ref #
25g (0.50) x 1-1/2" (38mm)	4004



RETROBULBAR ATKINSON ANESTHESIA NEEDLE

- Sterile, 10 per box


Description	Ref #
25g (0.50) x 1-1/2" (38mm)	4004B
23g (0.60) x 1-1/2" (38mm)	4004C



PERIBULBAR ATKINSON ANESTHESIA NEEDLE

- Sterile, 10 per box


Description	Ref #
25g (0.50) x 7/8" (22mm)	4004D
23g (0.60) x 1-1/4" (32mm)	4004F
25g (0.50) x 1-1/4" (32mm)	4004M



SUB TENON CANNULA

- Sterile, 10 per box

Description	Ref #
19g (1.10) x 1" (25mm) 25mm radius, flattened to 0.66mm	4005



SUB TENON CANNULA


- Sterile, 10 per box

Description	Ref #
19g (1.10)x 1" (25mm) Two 0.4mm opposite side ports	4005TP

CANNULAS, NEEDLES & CYSTOTOMES




IRRIGATION/AIR INJECTION & VISCOELASTIC




ANTERIOR CHAMBER CANNULA

- Sterile, 10 per box

Description	Ref #
25g (0.50) x 7/8" (22mm) 4mm bend	4001AF
27g (0.40) x 7/8" (22mm) 4mm bend	4002F
30g (0.30) x 7/8" (22mm) 4mm bend	4001F




Description	Ref #
30g (0.30) x 7/8" (22mm) 6mm bend	4001
25g (0.50) x 7/8" (22mm) 6mm bend	4001A
27g (0.40) x 7/8" (22mm) 6mm bend	4002




ANTERIOR CHAMBER/VISCOELASTIC CANNULA

- Sterile, 10 per box

Description	Ref #
23g (0.60) x 7/8" (22mm) 8mm bend, Thin Wall	4001CTW
27g (0.40) x 7/8" (22mm) 8mm bend	4002B



Description	Ref #
25g (0.50) x 7/8" (22mm) 10mm bend	4001AB
23g (0.60) x 7/8" (22mm) 10mm bend	4001B
27g (0.40) x 7/8" (22mm) 11mm bend	4002A




VISCOELASTIC CANNULA

- Sterile, 10 per box

Description	Ref #
27g (0.40) x 7/8" (22mm) 9mm bend, 45° Angle	4002L




IRRIGATION/AIR INJECTION & VISCOELASTIC



CURVED CANNULA

- Sterile, 10 per box

Description	Ref #
25g (0.50) x 1-1/8" (29mm)	4003




ANTERIOR CHAMBER WASHOUT CANNULA

- Sterile, 10 per box

Description	Ref #
19g (1.10) x 1" (25mm) 8mm bend	4036


CAPSULOTOMY




FORMED IRRIGATING CYSTOTOME

- Sterile, 10 per box

Description	Ref #
25g (0.50) x 7/8" (22mm)	4017
27g (0.40) x 7/8" (22mm)	4017B




Description	Ref #
25g (0.50) x 5/8" (16mm)	4017G
27g (0.40) x 5/8" (16mm)	4017IO
30g (0.30) x 5/8" (16mm)	4017I
25g (0.50) x 5/8" (16mm) 65° Angle	4017IE



STRAIGHT IRRIGATING CYSTOTOME

- Sterile, 10 per box


Description	Ref #
25g (0.50) x 7/8" (22mm)	4017A
27g (0.40) x 5/8" (16mm)	4017IS



REVERSE FORMED IRRIGATING CYSTOTOME

- Sterile, 10 per box


Description	Ref #
25g (0.50) x 5/8" (16mm)	4017L



ANGLED IRRIGATING CYSTOTOME

- Sterile, 10 per box


Description	Ref #
25g (0.50) x 5/8" (16mm) 12mm bend	4017M
27g (0.40) x 5/8" (16mm) 12mm bend	4017N



FORMED CAPSULORHEXIS CYSTOTOME

- Sterile, 10 per box
- Designed with cutting edge from front to back (push-pull) for circular tear. Tip is angled 90° and then rotated 90°

Description	Ref #
25g (0.40) x 5/8" (16mm)	4017S
27g (0.40) x 5/8" (16mm)	4017R




F-SL CURVE IRRIGATING CYSTOTOME

- Sterile, 10 per box

Description	Ref #
25g (0.40) x 5/8" (16mm) F-SL curve	4017Z
27g (0.40) x 5/8" F-SL curve	4017Y




HYDRODISSECTION



U-SHAPED HYDRODISSECTION CANNULA

- Sterile, 10 per box


Description	Ref #
25g (0.50) x 7/8" (22mm) 2.5mm L x 2.0mm W hook	4016AH



HORIZONTAL-FLOW HYDRODISSECTION CANNULA

- Sterile, 10 per box


Description	Ref #
25g (0.50) x 7/8" (22mm) 8mm bend, Flattened tip	4036A
27g (0.40) x 7/8" (22mm) 8mm bend, Flattened tip	4036B



VERTICAL-FLOW HYDRODISSECTION CANNULA

- Sterile, 10 per box


Description	Ref #
25g (0.50) x 7/8" (22mm) 4mm bend, Flattened tip	4036C
27g (0.40) x 7/8" (22mm) 4mm bend, Flattened tip	4036D



HELSINKI HYDRODISSECTION CANNULA

- Sterile, 10 per box
- Curved with flattened tip. (Horizontal Flow)


Description	Ref #
25g (0.50) x 7/8" (22mm) Curved, Flattened tip	4036F
27g (0.40) x 7/8" (22mm) Curved, Flattened tip	4036G



CHANG HYDRODISSECTION CANNULA

- Sterile, 10 per box
- 90° angled flattened tip provides a broad stream of fluid
- The point of the beveled tip is used to engage and rotate the nucleus to ensure complete separation


Description	Ref #
27g (0.40) x 3/4" (19mm) 1mm bend	4036J



SAUTER HYDRODISSECTION CANNULA

- Sterile, 10 per box
- Curved with flattened tip (Vertical Flow)


Description	Ref #
25g (0.50) x 7/8" (22mm)	4036M
27g (0.40) x 7/8" (22mm)	4036N




CORTICAL CLEAVING HYDRODISSECTOR

- Sterile, 10 per box

Description	Ref #
27g (0.40) x 7/8" (22mm) 11mm bend, round end	4036P



Description	Ref #
27g (0.40) x 7/8" (22mm) 11mm bend, flat end	4036Q




KELLAN HYDRODISSECTION CANNULA

- Sterile, 10 per box
- Flattened tip
- Curved, 20mm radius

Description	Ref #
25g (0.50) x 7/8" (22mm)	4068



CORTEX EXTRACTORS




SIMCOE CORTEX EXTRACTOR

- Sterile, 10 per box
- Curved, 0.3mm port

Description	Ref #
23g (0.60) x 1/2" (13mm)	4008

U-SHAPED CORTEX EXTRACTOR


- Sterile, 10 per box



Description	Ref #
21g (0.80) x 7/8" (22mm). Hook: 2.0mm L x 2.25mm W	4016
23g (0.60) x 7/8" (22mm). Hook: 2.0mm L x 2.25mm W	4016F
27g (0.40) x 7/8" (22mm). Hook: 2.5mm L x 2.0mm W	4016B

MICRO J-SHAPED CORTEX EXTRACTOR


- Sterile, 10 per box



Description	Ref #
25g (0.50) x 7/8" (22mm). Hook: 1.5mm L x 1.75mm W	4016G


CHARLEUX CORTEX EXTRACTOR

- Sterile, 10 per box
- Curved, 10mm radius, 45° tip



Description	Ref #
23g (0.60) x 1/2" (13mm)	4026


CAPSULE POLISHERS



KRATZ CAPSULE POLISHER

- Sterile, 10 per box
- Sand blasted tip


Description	Ref #
21g (0.80) x 7/8" (22mm) 9mm bend, 0.3mm top port	4027
23g (0.60) x 7/8" (22mm) 9mm bend, 0.3mm top port	4027A
21g (0.80) x 7/8" (22mm) 6mm bend, open end	4027B



OLIVE-TIP ANGLED CAPSULE POLISHER

- Sterile, 10 per box
- 2mm olive tip. Sand blasted tip


Description	Ref #
25g (0.50) x 1" (25mm) 10mm bend	4027E




TERRY CAPSULE POLISHER

- Sterile, 10 per box
- 10mm silicone sleeve

Description	Ref #
27g (0.40) x 7/8" (22mm) 6mm bend	4028
23g (0.60) x 7/8" (22mm) 7mm bend	4028B



IRRIGATING & ASPIRATING




SIMCOE DOUBLE I & A CANNULA

- Sterile, 10 per box.

Description	Ref #
19g/23g (0.60mm/0.60mm). Simcoe double I & A cannula, 0.30mm aspiration port	4007

LENS LOOP




LENS LOOP IRRIGATING VECTIS

- Sterile, 10 per box.
- For use with gravity fluid source.

Description	Ref #
24g (0.55) x 5.5mm (Three 0.25mm ports)	4029A

INFUSION




ANTERIOR CHAMBER MAINTAINER

- Sterile, 1 per pouch.
- 6" silicone tubing with 3mm tip
- Provides irrigation for chamber maintenance during anterior surgery

Description	Gauge	Ref #
Anterior Chamber Maintainer	20g (0.90)	4300P

INFUSION CANNULA


- Sterile, 10 per box.



Description	Ref #
2.0mm tip 20g x 2mm Long	4042-2S
4.0mm tip 20g x 4mm Long	4042-4S




LACRIMAL



CURVED OPEN END LACRIMAL CANNULA

- Sterile, 10 per box.


Description	Ref #
21g (0.80) x 1-1/2" (38mm)	6005



STRAIGHT CLOSED END LACRIMAL CANNULA

- Sterile, 10 per box.


Description	Ref #
21g (0.80) x 1-1/2" (38mm) 0.3mm port 1mm from tip	6006



CURVED CLOSED END LACRIMAL CANNULA

- Sterile, 10 per box.


Description	Ref #
25g (0.50) x 1-1/2" (38mm) 0.3mm port facing down 1mm from tip	6007



STRAIGHT STEPPED LACRIMAL CANNULA

- Sterile, 10 per box.

Description	Ref #
19g (1.10) x 1-1/2" (38mm) 12.5mm tip stepped down to 23g (0.60)	6008




CURVED OPEN END LACRIMAL CANNULA

- Sterile, 10 per box.

Description	Ref #
26g x 1-1/8", Curved, Open End	6010




LASIK CANNULAS



LASIK OPEN END CANNULA

- Sterile, 10 per box


Description	Ref #
25g (0.50) x 7/8" (22mm), formed 7mm from open end six 0.15mm side ports at 2mm, 4mm & 6mm from end	4225



LIN LASIK CANNULA

- Sterile, 10 per box


Description	Ref #
25g (0.50) x 7/8" (22mm), 9mm flattened angle, two 0.3mm ports	4036L



BURATTO LASIK CANNULA

- Sterile, 10 per box


Description	Ref #
25g (0.50) x 7/8" (22mm), 11mm angle, one top and two side ports	4036LB



LINDSTROM CANNULA

- Sterile, 10 per box

Description	Ref #
25g (0.50) x 7/8" (22mm) 45° angle, 5/8" radius, three 0.3mm ports	4036LR



SLADE LASIK IRRIGATOR

- Sterile, 10 per box

Description	Ref #
26g (0.45) x 1" (25mm) 9mm formed extension with spatulated tip	4036LS

SOFT CELL® PVA FOAM

*CE Marked



SOFT CELL PVA FOAM SURGICAL SPEAR

- Sterile, 150 per box.

Description	Ref #
Fiber-free, fast wicking. Firm	0525



SOFT CELL LASIK EYE DRAIN

- Sterile, 20 per box

Description	Ref #
(Chayet) 10mm internal diameter. Reduces tear film debris in the interface. Reduces reflux pooling at the hinge during ablation.	0522-5
(Chayet/Katzen) 11.5mm internal diameter. Reduces tear film debris in the interface. Reduces reflux pooling at the hinge during ablation.	0523-5



SOFT CELL EYE PLEDGET

- 150 per box

Description	Ref #
Eye Pledget	0527



SOFT CELL LASIK SHIELD

- 20 per box

Description	Ref #
9mm	0528




SOFT CELL WIPE + WICK

- 1 sterile wipe + wick per package, 20 packages per box

Description	Ref #
83mm x 83mm (Hydrated). Designed for use in delicate surgical procedures	0529



PUPIL EXPANSION DEVICES



*CE Marked



OASIS® IRIS EXPANDER

EXPAND YOUR VIEW. SUSTAIN YOUR ACCESS.

The OASIS® Iris Expander is a surgical solution to gently expand insufficiently dilated pupils during ophthalmic surgery and sustain visibility throughout the procedure.

The disposable device is indicated in cases where Miosis and/or Intraoperative Floppy Iris Syndrome are present.

Each OASIS® Iris Expander is a polypropylene ring packaged sterile with a loaded disposable inserter. Placement and removal are both gentle to fragile iris tissue.

PRODUCT FEATURES

- 6.25mm or 7.0mm diameters available
- Molded polypropylene ring
- Four pockets cradle delicate ocular tissue
- Packaged sterile with inserter
- Easy removal without potential tearing
- Single use

Product Name	Diameter	Ref #
OASIS Iris Expander 5pk	7.0mm	9700
OASIS Iris Expander 1pk	7.0mm	9700-S
OASIS Iris Expander 5pk	6.25mm	9625
OASIS Iris Expander 1pk	6.25mm	9625-S




OASIS® IRIS RETRACTORS

PRODUCT FEATURES

- Mechanical dilation of the pupil.
- Designed for ease of application.
- Smooth finish.
- Sterile, single-use.

AVAILABLE NOW

Product Name	Qty	Ref #
Iris Retractors, single use	30 per box	4465



*CE Marked

"In cases where the nucleus is 3 mm or larger and the pupil is less than 4 mm or under 5 mm after using an intraocular dilating agent, I plan from the outset to use the Iris Expander. This is true even in simple cases of small pupils or patients who won't dilate at the time of surgery. I find that most patients lose between 1 and 2 mm of dilation during surgery, so a 5-mm pupil can end up at 3 mm or less during cortical removal and lens insertion. The Iris Expander reduces the complication rate that accompanies cataract surgery in these small pupil cases."

Jeffery Whitman, MD, Key-Whitman Eye Center



ORDERING INFORMATION

Orders can be placed 24 hours via our website, fax, or the phone message center.

CUSTOMER SERVICE HOURS

7:00 AM - 4:30 PM (PST),

Monday - Friday, Excluding holidays.

USA: 844.820.8940

800.631.7210 FAX

INTERNATIONAL: 909.305.5400

909.305.9987 FAX

international@oasismedical.com

PAYMENT TERMS: NET 30 DAYS

SHIPPING: FOB GLENDORA, CALIFORNIA

Insurance and shipping charges added to invoice.

WEB (US ONLY): www.oasismedical.com

RETURN MERCHANDISE AUTHORIZATION (RMA)

Unopened items returned for credit must be accompanied by an RGA (return goods authorization number) and the reason for the return. An RMA can be requested through customer service.

Allowance Schedule for acceptable returns	
Sellable merchandise returned in factory condition within 60 days	Full Credit
Sellable merchandise returned in factory condition within 60-90 days	Cost less 40%
Sellable merchandise returned in factory condition 90 days to 1 year	Cost less 60%
Merchandise returned one year after invoice date	No Credit
Merchandise unsuitable for resale	No Credit
Merchandise defective at time of receipt by customer	Full Credit
Merchandise shipped in error	Full Credit

Merchandise that will not be accepted:

- “Special Order” or “Custom Products”
- Merchandise not in standard OASIS® packaging
- Opened products
- Discontinued products

NOTE: Customer must prepay all freight charges on stock returns except for defective material or merchandise shipped in error.

DOMESTIC SHIPPING POLICY

Shipments from OASIS® are FedEx only. Choices are Ground, 3-Day Select, 2-Day, Priority Overnight (by 10:30am) or Standard Overnight (by 4:30pm). Shipping costs will be added when order is placed. Customer’s FedEx or UPS account number may be used if provided to Customer Service at time order is placed. We cannot ship to PO boxes. We cannot be responsible for delays in transit. In-stock product orders received before 2pm PST Monday-Friday, excluding holidays, will usually ship on the same business day. Prepayment orders will ship once payment is received. Custom orders will require additional shipment times and will ship when the product is available.

WEBSITE ORDERS

Visit www.oasismedical.com to see terms of use, privacy policy, and returns authorization by web.



QUALITY ASSURANCE WORLDWIDE

The OASIS® Medical Quality System is structured to meet the requirements of the Food and Drug Administration (FDA) Quality System Regulations. OASIS® Medical is registered with the FDA as a:

- Medical Device Establishment - OASIS® Quality System complies with the Code of Federal Regulations (21 CFR §820)

The OASIS® Medical Quality Management System is certified by DQS Medizinprodukte GmbH to EN ISO 13485:2016, ISO 13485:2016, Medical Device Single Audit Program (MDSAP) and meets the requirements of Annex II - Excluding Section 4, of Council Directive 93/42/EEC concerning medical devices for applying the CE marking to products:



CE Marked Product Lines

Silicone Punctum Plugs	Extended Duration Plugs
Disposable Microsurgical Instruments	Surgical Marking Pen
HPMC Viscoelastic (International)	Microsurgical Knives
PVA Foam Sponge Products	Collagen Shields
Hydrogel Plugs	Iris Expander
Lacrimal Instruments	Iris Retractor

- Drug Establishment - OASIS® Quality System complies with the Code of Federal Regulations (21 CFR §211)

STERILIZATION METHOD

Products are sterilized by ethylene oxide, gamma radiation, electron beam, or steam according to EN 556 (European) guidelines and ANSI/AAMI/ISO 11135, 11137 or 17665 standards for medical device sterilization as applicable, to a sterility assurance level (SAL) of 10⁻⁶.

PRODUCT LABELING

OASIS® uses internationally recognized symbols on its product labeling.

	European Union CE Mark and notified body ID number		Manufacturer
	Consult instructions for use.		Do not resterilize
Rx Label	Physician use only		Do not use if package is damaged
	Do not reuse, single use		Keep dry
LS0812A	Lot or batch code		Temperature limit
	Serial Number		
yyyy-mm-dd	Use by. Placed next to Expiration Date – year/month		
[Part No.]	Reference Catalog Number (Numbers/Letter represent part number i.e. 6612D)		
	European Authorized Representative		

METHODS OF STERILIZATION

	Gamma or E-beam Radiation
	Ethylene Oxide
	Steam or Dry Heat



"Oasis TEARS became the dry eye drop of choice at our practice for a simple reason: our patients prefer them. Years ago, our team of physicians handed out many different artificial tears samples and asked patients for their feedback. Consistently, they gave the highest marks to Oasis TEARS.

Years later, Oasis TEARS are part of a comprehensive Premium Eye Care Package offered to our dry eye patients. Patients have been happy with the results and take comfort in knowing this is a high-quality product available only from their doctor's office."

Edward Lu, Director of Clinical Operations, IQ Laser Vision

Learn more and ask your OASIS® representative about our ocular surface solutions.




Donawa Lifescience Consulting Srl
Piazza Albania 10
00153 Rome, Italy



Glendora, CA, USA

(844) 820-8940 (USA Toll Free) | (909) 305-5400

customerservice@oasismedical.com | www.oasismedical.com



LIT-SRG-CAT Rev3 11-2019