



OASIS | SURGICAL®

PRODUCT CATALOG



844.820.8940

customerservice@oasismedical.com

oasismedical.com



OASIS | SURGICAL®



 **SOFT SHIELD®**
COLLAGEN SHIELD

CORNEAL COMFORT & PROTECTION

Provides ocular surface protection following surgery, injection and traumatic or non-traumatic corneal conditions. Available in Quick Shield, 12 hour, 24 hour, and 72 hour dissolution times.

Diameter: 14.05mm
Compound base curve: 9mm

PRODUCT FEATURES

- Comfort of a thin and pliable collagen barrier
- Collagen certified suitable for use by EDQM
- Hydrates to 70% water content
- Protects the surface of the eye
- Lubricates with collagen while it gradually dissolves
- CE approved

*CE Marked

Ref #	Product Name	Duration
7002	SOFT SHIELD Quick Shield Collagen Shield 10pk	Less than 12 hours
7012	SOFT SHIELD Collagen Shield 10pk	12 hours non-patched, 24 hours patched
7024	SOFT SHIELD Collagen Shield 10pk	24 hours non-patched, 72 hours patched
7072	SOFT SHIELD Collagen Shield 10pk	72 hours non-patched



 **VISCO SHIELD®**
TOPICAL DROPS

VISCO SHIELD® Topical Drops is preservative-free and designed to keep the ocular surface lubricated and hydrated to prevent drying of the eye.

PRODUCT FEATURES

- 2ml filled syringe with cannula
- Hypromellose
- Stores at room temperature

Ref #	Description	QTY
TD5200-0	5,000 cps	6 per box

SURGICAL CUSTOMIZED SOLUTIONS

SURGEON'S CHOICE CARD

Surgical Team - one peel pack reduces repetitive motion and prep time pulling disposable instruments

Surgery Center - supports efficiency with a continuous surgical flow in the O.R.

Surgeons - their preferred disposable instruments custom organized and on-hand


OASIS® Surgeon's Choice Card and Surgeon's Choice Cannula Card are alternative solutions to standard surgical preference cards. They provide the doctor's preferred surgical disposable instruments; such as any combination of up to eight OASIS cannulas, needles, and cystotomes (CNC), and up to two safety blades in an organized and compact delivery method. The Surgeon's Choice Cannula Card holds up to six CNC.

Both Choice Cards are packaged in a single peel pouch to reduce repetitive motion.

Surgical teams are able to reduce preparation time before surgery, reduce waste, and support a continuous flow from surgery to surgery

Contact your local OASIS Sales Consultant to request a quote or email: customsurgicalkit@oasismedical.com or contact Customer Service at **(844) 820-8940**

OPTION A



OPTION B





BULK PACKAGING

- Economical – reduced price per unit
- Saves time – no need to open boxes
- Less storage space needed
- Less waste disposal
- Provided sterile
- Automatic monthly shipments available

Available for cannulas, needles, cystotomes and Premier Edge® scalpels.

Surgical disposable products are now offered in this cost-effective and convenient Customized Inventory Management Program. The product is provided bulk and sterile in plastic bags, to eliminate waste and save time. Automatic monthly shipments are available with your choice of delivery each month.

Contact your local OASIS Sales Consultant to request a quote or email: customsurgicalkit@oasismedical.com

SCALPELS & KNIVES

OASIS® carries a wide variety of high quality ophthalmic microsurgical knives from Feather® and Premier Edge® for use in cataract, refractive, retinal and glaucoma surgery.



Feather redefines sharp and combines it with elegantly designed precision in their broad line of scalpels. Feather reusable scalpels are designed for repeated use, can be resterilized at your surgical facility in compliance with CMS regulations, and are easily disposed of after multiple uses. Feather also offers both single use disposable and safety scalpels.

Reusable	R	5 per box	Sure-grip aluminum handle	Stainless steel micro blades	Sterile individual autoclavable cases	
Disposable	D	5 per box	Disposable plastic handle	Stainless steel micro blades	Sterile blister pack	
Safeshield	S	5 per box	Disposable plastic handle	Stainless steel micro blades	Sterile blister pack	Retractable safety guard



***CE Marked**

Premier Edge Microsurgical Knives provide trustworthy performance during surgical procedures with enhanced precision on the cutting tip and ergonomic handle design. Safety scalpels are designed for one-handed use.

Disposable	D	6 per box	Disposable plastic handle	Stainless steel micro blades	Sterile blister pack	
Safeshield	S	6 per box	Disposable plastic handle	Stainless steel micro blades	Sterile blister pack	Retractable safety guard

GROOVE

R



FEATHER CONTROLLED DEPTH GROOVE

FEATHER® MICRO FEATHER Ophthalmic Scalpel with Aluminum Handle

- Designed for controlled depth incisions
- Sterile, 5 per box

Reusable	Depth
07330G	300 microns (0.3mm)
07335G	350 microns (0.35mm)
07340G	400 microns (0.4mm)
07355G	550 microns (0.55mm)
07360G	600 microns (0.6mm)

R



FEATHER GLAUCOMA GROOVE

FEATHER® MICRO FEATHER Ophthalmic Scalpel with Aluminum Handle

- Designed for fixed cutting depth
- Sterile, 5 per box

Reusable	Depth	Thickness
06340G	400 microns (0.4mm)	0.1mm

IMPLANT

R



FEATHER WOUND ENLARGEMENT IMPLANT

FEATHER® MICRO FEATHER Ophthalmic Scalpel with Aluminum Handle

- Sterile, 5 per box

Reusable	Width	Bevel	Blade
07635BR	3.5mm	Bevel Up	Angled
07651BR	5.1mm	Bevel Up	Angled

D



PREMIER EDGE IMPLANT

- Sterile, 6 per box

Disposable	Width	Bevel	Blade
PE3440	4.0mm	Bevel Up	Angled

SLIT KERATOME

R



SECTION VIEW

FEATHER SLIT KNIFE RIB & R

FEATHER® MICRO FEATHER Ophthalmic Scalpel with Aluminum Handle

- Sterile, 5 per box
- Improved Operability: Length from blade tip to bend is shortened from 8mm (conventional model) to 6mm.
- Arrowhead-shaped curved edge reduces incision resistance.
- Both sides of knife are thinly grinded (RIB shape).

Reusable	Disposable	Width	Bevel
0610RC		1.0mm	Bevel Up
0611RC		1.1mm	Bevel Up
0612RC		1.2mm	Bevel Up
0614RC	P-0614RC	1.4mm	Bevel Up
0616RC	P-0616RC	1.6mm	Bevel Up
0618RC	P-0618RC	1.8mm	Bevel Up
0620RC	P-0620RC	2.0mm	Bevel Up
0622RC	P-0622RC	2.2mm	Bevel Up
0624RC	P-0624RC	2.4mm	Bevel Up
06275RC	P-0627RC	2.7mm	Bevel Up
0628RC	P-0628RC	2.8mm	Bevel Up
0630RC	P-0630RC	3.0mm	Bevel Up

R D



FEATHER ULTRA-SHARP SLIT/KERATOME

FEATHER® MICRO FEATHER Ophthalmic Scalpel with Aluminum Handle

FEATHER® MICRO FEATHER Disposable Ophthalmic Scalpel

Reusable	Disposable	Width	Bevel	Blade
0614BUC	P-0614BUC	1.4mm	Bevel Up	Angled
0616BUC	P-0616BUC	1.6mm	Bevel Up	Angled
0618BUC	P-0618BUC	1.8mm	Bevel Up	Angled
0622BUC	P-0622BUC	2.2mm	Bevel Up	Angled
0624BUC	P-0624BUC	2.4mm	Bevel Up	Angled
0625BUC		2.5mm	Bevel Up	Angled
06265BUC		2.6mm	Bevel Up	Angled
06275BUC	P-06275BUC	2.75mm	Bevel Up	Angled
0628BUC	P-0628BUC	2.8mm	Bevel Up	Angled
0630BUC	P-0630BUC	3.0mm	Bevel Up	Angled
0632BUC	P-0632BUC	3.2mm	Bevel Up	Angled

D



PREMIER EDGE SLIT

• Sterile, 6 per box

Disposable	Width	Bevel	Blade
PE3818	1.8mm	Bevel Up	Angled
PE3822	2.2mm	Bevel Up	Angled
PE3824	2.4mm	Bevel Up	Angled
PE3825	2.5mm	Bevel Up	Angled
PE3826	2.65mm	Bevel Up	Angled
PE3827	2.75mm	Bevel Up	Angled
PE3828	2.85mm	Bevel Up	Angled
PE3830	3.0mm	Bevel Up	Angled
PE3832	3.2mm	Bevel Up	Angled
PE3832S	3.2mm	Bevel Up	Angled

SLIT KERATOME

D




PREMIER EDGE THIN SLIT

• Sterile, 6 per box

Disposable	Safety	Width	Bevel	Blade
PE4820-TL	PE4820-SL	2.0mm	Bevel Up	Angled
PE4822-TL	PE4822-SL	2.2mm	Bevel Up	Angled
PE4824-TL	PE4824-SL	2.4mm	Bevel Up	Angled
PE4826-TL	PE4826-SL	2.65mm	Bevel Up	Angled
PE4827-TL	PE4827-SL	2.75mm	Bevel Up	Angled
PE4828-TL	PE4828-SL	2.85mm	Bevel Up	Angled
PE4830-TL	PE4830-SL	3.0mm	Bevel Up	Angled
PE4832-TL	PE4832-SL	3.2mm	Bevel Up	Angled

S



D



PREMIER EDGE CLEAR CORNEA SLIT

• Sterile, 6 per box

Disposable	Safety	Width	Bevel	Blade
PE3122		2.0mm	Double	Angled
PE3125		2.5mm	Double	Angled
PE3126	PE3126-SL	2.65mm	Double	Angled
PE3127	PE3127-SL	2.75mm	Double	Angled
PE3128	PE3128-SL	2.85mm	Double	Angled
PE3132		3.2mm	Double	Angled

D



PREMIER EDGE CLEAR CORNEA PARALLEL SLIT

• Sterile, 6 per box

Disposable	Width	Bevel	Blade
PE3226	2.65mm	Double	Angled
PE3227	2.75mm	Double	Angled
PE3228	2.85mm	Double	Angled
PE3230	3.0mm	Double	Angled

D



PREMIER EDGE CLEAR CORNEA TRAPEZOID SLIT

• Sterile, 6 per box

Disposable	Width Mark	Bevel	Blade
PE3303	1.5mm - 1.7mm	Bevel Up	Angled
PE3301	1.8mm - 2.2mm	Bevel Up	Angled
PE3302	2.5mm - 3.0mm	Bevel Up	Angled

MINI BLADES

D



PREMIER EDGE MINI BLADES

• Sterile

Disposable	Blade
PE5700	Multi-sided, sharp all around, 6 per box
PE6400	Rounded Tip, sharp on one side, 20 per box
PE6700	Curved Tip, sharp on one side, 20 per box
PE6900	Rounded Tip, sharp all around, 6 per box

CRESCENT TUNNEL

R D



FEATHER CRESCENT TUNNEL

FEATHER® MICRO FEATHER Ophthalmic Scalpel with Aluminum Handle

FEATHER® MICRO FEATHER Disposable Ophthalmic Scalpel

- 45° angle
- Sterile, 5 per box

Reusable	Disposable	Width	Bevel	Blade
07922B		2.2mm	Double	Angled
07926BU	P-7926BU	2.6mm	Bevel Up	Angled
07926BD		2.6mm	Bevel Down	Angled
07926B	P-7926B	2.6mm	Double	Angled

D S



PREMIER EDGE CRESCENT TUNNEL

- Sterile, 6 per box

Disposable	Safety	Width	Bevel	Blade
PE3720A	PE3720A-SL	2.0mm	Bevel Up	Angled
PE3720		2.0mm	Bevel Up	Straight

SCLERAL

R



FEATHER MULTISIDED SCLERAL

FEATHER® MICRO FEATHER Ophthalmic Scalpel with Aluminum Handle

- Sterile, 5 per box

Reusable	Bevel	Thickness
0780	Double	0.38mm

ROUND TUNNEL

R D S



FEATHER ROUND TUNNEL

FEATHER® MICRO FEATHER Ophthalmic Scalpel with Aluminum Handle
 FEATHER® MICRO FEATHER Disposable Ophthalmic Scalpel
 FEATHER® MICRO FEATHER SAFESHIELD Disposable Ophthalmic Scalpel
 • Sterile, 5 per box • 0.25mm blade thickness

Reusable	Disposable	Width	Bevel	Blade	Thickness
07722B		2.2mm	Double	Angled	0.25mm
07726B	P-7726B	2.6mm	Double	Angled	0.25mm
07730B		3.0mm	Double	Angled	0.25mm

D



PREMIER EDGE ROUND TUNNEL

• Sterile, 6 per box

Disposable	Width	Bevel	Blade
PE3622	2.2mm	Double	Angled
PE3630	3.0mm	Double	Angled

"The delivery of first-rate surgical outcomes is of paramount importance to our team. We constantly strive to produce those outcomes in the most efficient manner possible. Cost containment in our supply chain allows us to conserve our resources and to redirect them where they matter most. This is critical for any ambulatory surgery center to thrive. We use FEATHER Reusable Scalpels as one means to achieve these goals. The product is superiorly made and this allows us, along with the proper care, to get the most use from each blade."

Anonymous Independent ASC

STAB INCISION

R D S



FEATHER STAB INCISION

FEATHER® MICRO FEATHER Ophthalmic Scalpel with Aluminum Handle

FEATHER® MICRO FEATHER Disposable Ophthalmic Scalpel

FEATHER® MICRO FEATHER SAFESHIELD Disposable Ophthalmic Scalpel

• Sterile, 5 per box

Reusable	Disposable	Safety	Angle	Blade
0715	P-715	PSS-715	15°	Straight
0722	P-722		22.5°	Straight
0730	P-730	PSS-730	30°	Straight
0745	P-745		45°	Straight

D



PREMIER EDGE STAB INCISION

• Sterile, 6 per box

Disposable	Safety	Angle	Bevel	Blade
PE3015	PE3015-SL	15°	Double	Straight
PE3022		22.5°	Double	Straight
PE3030		30°	Double	Straight
PE3045		45°	Double	Straight

S



PREMIER EDGE RESTRICTED DEPTH STAB INCISION

• Sterile, 6 per box

Disposable	Angle	Depth
PE3015-3	15°	3.0mm
PE3015-5	15°	5.0mm

D



MVR

R



FEATHER MVR

FEATHER® MICRO FEATHER Ophthalmic Scalpel with Aluminum Handle

- Sterile, 5 per box

Reusable	Disposable	Gauge	Bevel	Blade
07520	P-7520B	20g	Double	Straight
07523		20g	Double	Straight

D



PREMIER EDGE MVR

- Sterile, 6 per box

Disposable	Safety	Gauge	Bevel	Blade
PE3519		19g	Double	Straight
PE3520	PE3520-SL	20g	Double	Straight
PE3523		23g	Double	Straight
PE3525		25g	Double	Straight
PE3519A		19g	Double	Angled
PE3520A		20g	Double	Angled
PE3523A		23g	Double	Angled

S



SIDE PORT

D



PREMIER EDGE SIDE PORT

- Sterile, 6 per box

Disposable	Width	Bevel	Blade
PE3010	1.0mm	Double	Angled
PE3012	1.2mm	Double	Angled
PE3014	1.4mm	Double	Angled

R



FEATHER SIDE PORT

FEATHER® MICRO FEATHER Ophthalmic Scalpel with Aluminum Handle

- Sterile, 5 per box

Reusable	Width	Bevel
0610RC	1.0mm	Bevel Up
0611RC	1.1mm	Bevel Up
0612RC	1.2mm	Bevel Up


SECTION VIEW

DISPOSABLE MICROSURGICAL INSTRUMENTS

OASIS® supplies an extensive line of cannulas, needles, and cystotomes.

*CE Marked


ANESTHESIA



RETROBULBAR ANESTHESIA NEEDLE

- Sterile, 10 per box


Ref #	Description
4004	25g (0.50) x 1-1/2" (38mm) 4004



RETROBULBAR ATKINSON ANESTHESIA NEEDLE

- Sterile, 10 per box


Ref #	Description
4004B	25g (0.50) x 1-1/2" (38mm)
4004C	23g (0.60) x 1-1/2" (38mm)



PERIBULBAR ATKINSON ANESTHESIA NEEDLE

- Sterile, 10 per box


Ref #	Description
4004D	25g (0.50) x 7/8" (22mm)
4004F	23g (0.60) x 1-1/4" (32mm)
4004M	25g (0.50) x 1-1/4" (32mm)



SUB TENON CANNULA

- Sterile, 10 per box

Ref #	Description
4005	19g (1.10) x 1" (25mm) 25mm radius, flattened to 0.66mm




SUB TENON CANNULA

- Sterile, 10 per box

Ref #	Description
4005TP	19g (1.10) x 1" (25mm) Two 0.4mm opposite side ports


IRRIGATION/AIR INJECTION & VISCOELASTIC



ANTERIOR CHAMBER CANNULA

- Sterile, 10 per box


Ref #	Description
4001AF	25g (0.50) x 7/8" (22mm) 4mm bend
4002F	27g (0.40) x 7/8" (22mm) 4mm bend
4001F	30g (0.30) x 7/8" (22mm) 4mm bend



ANTERIOR CHAMBER CANNULA

- Sterile, 10 per box


Ref #	Description
4001	30g (0.30) x 7/8" (22mm) 6mm bend
4001A	25g (0.50) x 7/8" (22mm) 6mm bend
4002	27g (0.50) x 7/8" (22mm) 6mm bend



ANTERIOR CHAMBER/VISCOELASTIC CANNULA

- Sterile, 10 per box



Ref #	Description
4001CTW	23g (0.60) x 7/8" (22mm) 8mm bend, Thin Wall
4002B	27g (0.40) x 7/8" (22mm) 8mm bend



ANTERIOR CHAMBER/VISCOELASTIC CANNULA

- Sterile, 10 per box

Ref #	Description
4001AB	25g (0.50) x 7/8" (22mm) 10mm bend
4001B	23g (0.60) x 7/8" (22mm) 10mm bend
4002A	27g (0.40) x 7/8" (22mm) 11mm bend




VISCOELASTIC CANNULA

- Sterile, 10 per box

Ref #	Description
4002L	27g (0.40) x 7/8" (22mm) 9mm bend, 45° Angle


IRRIGATION/AIR INJECTION & VISCOELASTIC



CURVED CANNULA

- Sterile, 10 per box

Ref #	Description
4003	25g (0.50) x 1-1/8" (29mm)




ANTERIOR CHAMBER WASHOUT CANNULA

- Sterile, 10 per box

Ref #	Description
4036	19g (1.10) x 1" (25mm) 8mm bend


CAPSULOTOMY



FORMED IRRIGATING CYSTOTOME

- Sterile, 10 per box


Ref #	Description
4017	25g (0.50) x 7/8" (22mm)
4017B	27g (0.40) x 7/8" (22mm)



FORMED IRRIGATING CYSTOTOME

- Sterile, 10 per box


Ref #	Description
4017G	25g (0.50) x 5/8" (16mm)
4017IO	27g (0.40) x 5/8" (16mm)
4017I	30g (0.30) x 5/8" (16mm)
4017IE	25g (0.50) x 5/8" (16mm) 65° Angle



STRAIGHT IRRIGATING CYSTOTOME

- Sterile, 10 per box


Ref #	Description
4017A	25g (0.50) x 7/8" (22mm)
4017IS	27g (0.40) x 5/8" (16mm)



REVERSE FORMED IRRIGATING CYSTOTOME

- Sterile, 10 per box


Ref #	Description
4017L	25g (0.50) x 5/8" (16mm)



ANGLED IRRIGATING CYSTOTOME

- Sterile, 10 per box


Ref #	Description
4017M	25g (0.50) x 5/8" (16mm) 12mm bend
4017N	27g (0.40) x 5/8" (16mm) 12mm bend



FORMED CAPSULORHEXIS CYSTOTOME

- Sterile, 10 per box
- Designed with cutting edge from front to back (push-pull) for circular tear. Tip is angled 90° and then rotated 90°

Ref #	Description
4017S	25g (0.40) x 5/8" (16mm)
4017R	27g (0.40) x 5/8" (16mm)




F-SL CURVE IRRIGATING CYSTOTOME

- Sterile, 10 per box

Ref #	Description
4017Z	25g (0.40) x 5/8" (16mm) F-SL curve
4017Y	27g (0.40) x 5/8" F-SL curve


HYDRODISSECTION



U-SHAPED HYDRODISSECTION CANNULA

- Sterile, 10 per box


Ref #	Description
4016AH	25g (0.50) x 7/8" (22mm) 2.5mm L x 2.0mm W hook



HORIZONTAL-FLOW HYDRODISSECTION CANNULA

- Sterile, 10 per box


Ref #	Description
4036A	19g (1.10) x 1" (25mm) 8mm bend
4036B	27g (0.40) x 7/8" (22mm) 8mm bend, Flattened tip



VERTICAL-FLOW HYDRODISSECTION CANNULA

- Sterile, 10 per box


Ref #	Description
4036C	25g (0.50) x 7/8" (22mm) 4mm bend, Flattened tip
4036D	27g (0.40) x 7/8" (22mm) 4mm bend, Flattened tip



VERTICAL-FLOW HYDRODISSECTION CANNULA

- Sterile, 10 per box
- Curved with flattened tip. (Horizontal Flow)


Ref #	Description
4036F	25g (0.50) x 7/8" (22mm) Curved, Flattened tip
4036G	27g (0.40) x 7/8" (22mm) Curved, Flattened tip



CHANG HYDRODISSECTION CANNULA

- Sterile, 10 per box
- 90° angled flattened tip provides a broad stream of fluid
- The point of the beveled tip is used to engage and rotate the nucleus to ensure complete separation


Ref #	Description
4036J	27g (0.40) x 3/4" (19mm) 1mm bend



SAUTER HYDRODISSECTION CANNULA

- Sterile, 10 per box
- Curved with flattened tip (Vertical Flow)


Ref #	Description
4036M	25g (0.50) x 7/8" (22mm)
4036N	27g (0.40) x 7/8" (22mm) 4036N



CORTICAL CLEAVING HYDRODISSECTOR

- Sterile, 10 per box


Ref #	Description
4036P	27g (0.40) x 7/8" (22mm) 11mm bend, round end



CORTICAL CLEAVING HYDRODISSECTOR

- Sterile, 10 per box

Ref #	Description
4036Q	27g (0.40) x 7/8" (22mm) 11mm bend, flat end




KELLAN HYDRODISSECTION CANNULA

- Sterile, 10 per box
- Flattened tip
- Curved, 20mm radius

Ref #	Description
4068	25g (0.50) x 7/8" (22mm)


CORTEX EXTRACTORS



SIMCOE CORTEX EXTRACTOR

- Sterile, 10 per box
- Curved, 0.3mm port


Ref #	Description
4008	23g (0.60) x 1/2" (13mm)



U-SHAPED CORTEX EXTRACTOR

- Sterile, 10 per box


Ref #	Description
4016	21g (0.80) x 7/8" (22mm). Hook: 2.0mm L x 2.25mm W
4016F	23g (0.60) x 7/8" (22mm). Hook: 2.0mm L x 2.25mm W
4016B	27g (0.40) x 7/8" (22mm). Hook: 2.5mm L x 2.0mm W



MICRO J-SHAPED CORTEX EXTRACTOR

- Sterile, 10 per box

Ref #	Description
4016G	25g (0.50) x 7/8" (22mm). Hook: 1.5mm L x 1.75mm W





CHARLEUX CORTEX EXTRACTOR

- Sterile, 10 per box
- Curved, 10mm radius, 45° tip

Ref #	Description
4026	23g (0.60) x 1/2" (13mm)


CAPSULE POLISHERS



KRATZ CAPSULE POLISHER

- Sterile, 10 per box
- Sand blasted tip


Ref #	Description
4027	21g (0.80) x 7/8" (22mm) 9mm bend, 0.3mm top port
4027A	23g (0.60) x 7/8" (22mm) 9mm bend, 0.3mm top port
4027B	21g (0.80) x 7/8" (22mm) 6mm bend, open end



OLIVE-TIP ANGLED CAPSULE POLISHER

- Sterile, 10 per box
- 2mm olive tip. Sand blasted tip

Ref #	Description
4027E	25g (0.50) x 1" (25mm) 10mm bend




TERRY CAPSULE POLISHER

- Sterile, 10 per box
- 10mm silicone sleeve

Ref #	Description
4028	27g (0.40) x 7/8" (22mm) 6mm bend
4028B	23g (0.60) x 7/8" (22mm) 7mm bend

IRRIGATING & ASPIRATING




SIMCOE DOUBLE I & A CANNULA

- Sterile, 10 per box

Ref #	Description
4007	19g/23g (0.60mm/0.60mm). Simcoe double I & A cannula, 0.30mm aspiration port

LENS LOOP




LENS LOOP IRRIGATING VECTIS

- Sterile, 10 per box
- For use with gravity fluid source

Ref #	Description
4029A	24g (0.55) x 5.5mm (Three 0.25mm ports)


INFUSION



ANTERIOR CHAMBER MAINTAINER

- Sterile, 1 per pouch.
- 6" silicone tubing with 3mm tip
- Provides irrigation for chamber maintenance during anterior surgery

Ref #	Gauge	Description
4029A	20g (0.90)	Anterior Chamber Maintainer





INFUSION CANNULA

- Sterile, 10 per box

Ref #	Description
4042-2S	2.0mm tip 20g x 2mm Long
4042-4S	4.0mm tip 20g x 4mm Long


LACRIMAL



CURVED OPEN END LACRIMAL CANNULA

• Sterile, 10 per box


Ref #	Description
6005	21g (0.80) x 1-1/2" (38mm)



STRAIGHT CLOSED END LACRIMAL CANNULA

• Sterile, 10 per box


Ref #	Description
6006	21g (0.80) x 1-1/2" (38mm) 0.3mm port 1mm from tip



CURVED CLOSED END LACRIMAL CANNULA

• Sterile, 10 per box



Ref #	Description
6007	25g (0.50) x 1-1/2" (38mm) 0.3mm port facing down 1mm from tip



STRAIGHT STEPPED LACRIMAL CANNULA

• Sterile, 10 per box

Ref #	Description
6008	19g (1.10) x 1-1/2" (38mm) 12.5mm tip stepped down to 23g (0.60)




CURVED OPEN END LACRIMAL CANNULA

• Sterile, 10 per box

Ref #	Description
6010	26g x 1-1/8", Curved, Open End


LASIK CANNULAS



LASIK OPEN END CANNULA

- Sterile, 10 per box


Ref #	Description
4225	25g (0.50) x 7/8" (22mm), formed 7mm from open end six 0.15mm side ports at 2mm, 4mm & 6mm from end



LIN LASIK CANNULA

- Sterile, 10 per box


Ref #	Description
4036L	25g (0.50) x 7/8" (22mm), 9mm flattened angle, two 0.3mm ports



BURATTO LASIK CANNULA

- Sterile, 10 per box


Ref #	Description
4036LB	25g (0.50) x 7/8" (22mm), 11mm angle, one top and two side ports



LINDSTROM CANNULA

- Sterile, 10 per box

Ref #	Description
4036LR	25g (0.50) x 7/8" (22mm) 45° angle, 5/8" radius, three 0.3mm ports



SLADE LASIK IRRIGATOR

- Sterile, 10 per box

Ref #	Description
4036LS	26g (0.45) x 1" (25mm) 9mm formed extension with spatulated tip

SOFT CELL® PVA FOAM



SOFT CELL PVA FOAM SURGICAL SPEAR

- Sterile, 150 per box







Ref #	Description
0525	Fiber-free, fast wicking. Firm

*CE Marked

MARKETING MATERIALS


OASIS® offers educational materials in the form of brochures, displays, posters, and videos for use in your practice and to distribute to your patients.

To view and download digital versions of all resources, go to Practice & Patient Materials at oasismedical.com



Ref #	Product Name		Ref #	Product Name	
LIT-USDE-CAT	Dry Eye Catalog		LIT-OTRX-W	Oasis TEARS Doctor Recommendation Pad, with OASIS Contact Information	
6544C	Oasis TEARS Brochure		LIT-OTRX-NB	Oasis TEARS Doctor Recommendation Pad, without OASIS Contact Information	
LIT-PP-PB	Punctal Plug Brochure		LIT-DES-WC	Dry Eye Symptoms Wall Chart	
LIT-DEP-PB (English) LIT-DEP-PBS (Spanish) LIT-DEP-PBC (Trad. Chinese)	Dry Eye Products Patient Brochure		LIT-WCDE	Punctal Plug Wall Chart	
LIT-OT-POP	Oasis TEARS Point-of-Purchase Display		LIT-SRG-CAT	Surgical Catalog	



PUPIL EXPANSION DEVICES



*CE Marked



OASIS® IRIS EXPANDER

EXPAND YOUR VIEW. SUSTAIN YOUR ACCESS.

The OASIS® Iris Expander is a surgical solution to gently expand insufficiently dilated pupils during ophthalmic surgery and sustain visibility throughout the procedure.

The disposable device is indicated in cases where Miosis and/or Intraoperative Floppy Iris Syndrome are present.

Each OASIS® Iris Expander is a polypropylene ring packaged sterile with a loaded disposable inserter. Placement and removal are both gentle to fragile iris tissue.

PRODUCT FEATURES

- 6.25mm or 7.0mm diameters available
- Molded polypropylene ring
- Four pockets cradle delicate ocular tissue
- Packaged sterile with inserter
- Easy removal without potential tearing
- Single use

Ref #	Product Name	Diameter
9700	OASIS Iris Expander 5pk	7.0mm
9700-S	OASIS Iris Expander 1pk	7.0mm
9625	OASIS Iris Expander 5pk	6.25mm
9625-S	OASIS Iris Expander 1pk	6.25mm

OPTALMED® IRIS RETRACTORS

PRODUCT FEATURES

- Mechanical dilation of the pupil.
- Designed for ease of application.
- Smooth finish.
- Sterile, single-use.

Ref #	Product Name	Qty
OPTALMED-7002	Optalmed Iris Retractors, single use	30 per box

ACCESSORIES & DRY EYE DIAGNOSTICS



6648



4003-IK



6004-27

Ref #	Product Name	Description	Qty Per Box
6512	Tear Flo Test Strips	Sterile strip	100

Ref #	Product Name	Description	Qty Per Box
6625	SOFT PLUG Punctal Plug Gauge	0.5mm and 0.6mm tips	1
6627	SOFT PLUG Punctal Plug Gauge	0.7mm and 0.8mm tips	1
6630	SOFT PLUG Punctal Plug Gauge Set	Includes both instruments and sterilization tray	1 set

Ref #	Product Name	Description	Qty Per Box
6643	Straight Jewelers Forceps	Straight #3	1
6647	Curved Jewelers Forceps	Curved #7	1
6648	Grooved Forceps	For the insertion of Extended Duration & Collagen Plugs	1
6640	Castroviejo Lacrimal Dilator	Double End	1

Ref #	Product Name	Description	Qty Per Box
4003-IK	Lacrimal Irrigation Kit	One 3mL syringe and one 25g x 1 1/8" lacrimal cannula, curved, open end	10
6001	Lacrimal Intubation Set	23g	3
6004	Lacrimal Intubation Set with Retrieval Device	23g	3
6004-27	Lacrimal Intubation Kit	27g	3

FITOVER SUNGLASSES



- Comfortable fit
- Lightweight frame
- Polycarbonate lenses
- Provides 48% more protection from sunlight than ordinary sunglasses
- 100% UVA/UVB protection

Ref #	Description
OM8100	OASIS Cataract Glasses

ORDERING INFORMATION

Orders can be placed 24 hours via our website, fax, or the phone message center.

CUSTOMER SERVICE HOURS
7:00 AM - 4:30 PM (PST),
Monday - Friday, Excluding holidays.

USA: 844.820.8940
800.631.7210 FAX

PAYMENT TERMS: NET 30 DAYS
SHIPPING: FOB GLENDORA, CALIFORNIA
Insurance and shipping charges added to invoice.

INTERNATIONAL: 909.305.5400
909.305.9987 FAX
international@oasismedical.com

WEB (US ONLY): myoasis360.com

RETURN MERCHANDISE AUTHORIZATION (RMA)

Unopened items returned for credit must be accompanied by an RGA (return goods authorization number) and the reason for the return. An RMA can be requested through customer service.

Allowance Schedule for acceptable returns	
Sellable merchandise returned in factory condition within 60 days	Full Credit
Sellable merchandise returned in factory condition within 60-90 days	Cost less 40%
Sellable merchandise returned in factory condition 90 days to 1 year	Cost less 60%
Merchandise returned one year after invoice date	No Credit
Merchandise unsuitable for resale	No Credit
Merchandise defective at time of receipt by customer	Full Credit
Merchandise shipped in error	Full Credit

Merchandise that will not be accepted:

- “Special Order” or “Custom Products”
- Merchandise not in standard OASIS® packaging
- Opened products
- Discontinued products

NOTE: Customer must prepay all freight charges on stock returns except for defective material or merchandise shipped in error.

DOMESTIC SHIPPING POLICY

Shipments from OASIS® are FedEx only. Choices are Ground, 3-Day Select, 2-Day, Priority Overnight (by 10:30am) or Standard Overnight (by 4:30pm). Shipping costs will be added when order is placed. Customer’s FedEx or UPS account number may be used if provided to Customer Service at time order is placed. We cannot ship to PO boxes. We cannot be responsible for delays in transit. In-stock product orders received before 2pm PST Monday-Friday, excluding holidays, will usually ship on the same business day. Prepayment orders will ship once payment is received. Custom orders will require additional shipment times and will ship when the product is available.

WEBSITE ORDERS

Visit myoasis360.com to see terms of use, privacy policy, and returns authorization by web.

FOR GENERAL INFO AND CATALOG FOR DOWNLOAD:

Visit oasismedical.com

QUALITY ASSURANCE WORLDWIDE

The OASIS® Medical Quality System is structured to meet the requirements of the Food and Drug Administration (FDA) Quality System Regulations. OASIS® Medical is registered with the FDA as a:

- Medical Device Establishment - OASIS® Quality System complies with the Code of Federal Regulations (21 CFR §820)

The OASIS® Medical Quality Management System is certified by DQS Medizinprodukte GmbH to EN ISO 13485:2016, Medical Device Single Audit Program (MDSAP) and meets the requirements of Annex II - Excluding Section 4, of Council Directive 93/42/EEC and Directive 2026/607 concerning medical devices for applying the CE marking to products:



CE Marked Product Lines

- | | |
|--------------------------------------|-------------------------|
| Silicone Punctum Plugs | Extended Duration Plugs |
| Disposable Microsurgical Instruments | Microsurgical Knives |
| HPMC Viscoelastic (International) | Collagen Shields |
| PVA Foam Spears | Iris Expander |
| Hydrogel Plugs | |

- Drug Establishment - OASIS® Quality System complies with the Code of Federal Regulations (21 CFR §211)

STERILIZATION METHOD

Products are sterilized by ethylene oxide, gamma radiation, electron beam, or steam according to EN 556 (European) guidelines and ANSI/AAMI/ISO 11135, 11137 or 17665 standards for medical device sterilization as applicable, to a sterility assurance level (SAL) of 10⁻⁶.

METHODS OF STERILIZATION

- | | |
|--|--|
| | Gamma or E-beam Radiation, Single Sterile Barrier System |
| | Gamma or E-beam Radiation, Double Sterile Barrier System |
| | Gamma or E-beam Radiation, Single Sterile Barrier System with Protective Packaging Inside |
| | Gamma or E-beam Radiation, Single Sterile Barrier System with Protective Packaging Outside |
| | Ethylene Oxide, Single Sterile Barrier System |
| | Ethylene Oxide, Double Sterile Barrier System |
| | Steam or Dry Heat, Single Sterile Barrier System |

PRODUCT LABELING

OASIS® uses internationally recognized symbols on its product labeling

- | | | | |
|--|--|--|---|
| | European Union CE Mark and notified body ID number | | Manufacturer |
| | Consult instructions for use. | | UDI |
| | Physician use only | | Patient Information Website |
| | Do not reuse, single use | | Contains biological material of animal origin |
| | Lot or batch code | | Swiss Authorized Representative |
| | Medical Device | | Do not resterilize |
| | Use by. Placed next to Expiration Date – year/month | | Do not use if package is damaged |
| | Reference Catalog Number (Numbers/Letter represent part number i.e. 6612D) | | No Latex |
| | European Authorized Representative | | Keep dry |
| | | | Temperature limit |



EC	REP
----	-----

Donawa Lifescience Consulting Srl
Piazza Albania 10
00153 Rome, Italy

CH	REP
----	-----

MedEnvoy Switzerland
Gotthardstrasse 28 6302 Zug
Switzerland

LIT-SRG-CAT Rev4 04-2025

OASIS 
2035 E Financial Way
Glendora, CA 91741